

18 Litre



N-CLASS

TEMP AND NAME OF CYCLE		STERILIZABLE LOADS ¹	TECHNICAL SPECS OF THE CYCLES ²				
				STR	DRY	TOT	KG
121° N	Solids Open	Solid load	●	18:00	15:00	~48m	1.5
134° N	Solids Open		●	06:00	15:00	~36m	1.5
134° N	Prion	Solid load related to the Creutzfeldt-Jakob disease	●	20:00	15:00	~50m	1.5
121° N	Rapid	Solid load	●	18:00	05:00	~38m	3.0
134° N	Rapid		●	06:00	05:00	~28m	3.0

REFERRING TO MIDMARK DENTISURE N18

N-CLASS

TEMP AND NAME OF CYCLE		STERILIZABLE LOADS ¹	TECHNICAL SPECS OF THE CYCLES ²				
				STR	DRY	TOT	KG
121° N	Solids Open	Solid load	●	18:00	20:00	~58m	2.5
134° N	Solids Open		●	06:00	20:00	~46m	2.5
134° N	Prion	Solid load related to the Creutzfeldt-Jakob disease	●	20:00	20:00	~60m	2.5
121° N	Rapid	Solid load	●	18:00	10:00	~48m	4.5
134° N	Rapid		●	06:00	10:00	~38m	4.5

REFERRING TO MIDMARK DENTISURE N23

23 Litre