



MEDICAL INTERNATIONAL MODELS

627/224 Barrier-Free® Examination Chairs



Accessible design helps ensure equal care for all.

Accessibility in healthcare goes beyond the ability of patients to access the care they need. The physical accessibility of a facility and the equipment within it plays a significant role. The proper exam chair is important in providing a thorough examination. Barrier-Free examination chairs by Midmark are designed to make accessing the exam chair easier, even for those with mobility concerns or other physical limitations.

Midmark 627 Barrier-Free® Examination Chair



01 BARRIER-FREE LOW CHAIR HEIGHT
can make accessing the examination chair easier for patients. (46 cm or 18 in minimum height)

02 POWERED BACK + HEIGHT
assists in positioning the patient for more comfortable examinations.

03 CLEAN ASSIST ROLLER SYSTEM
is a retractable roller base that allows you to easily and safely move the chair for cleaning.

04 STANDARD STIRRUPS
are integrated into the exam chair for lower-body examinations.

05 PELVIC TILT
makes proper exam positioning easier.

06 FRONT STORAGE DRAWER
gives you room to store supplies at the point of care.

07 ACTIVE SENSING TECHNOLOGY®
pauses chair movement when disrupting impact is detected.



Lithotomy Position
Easily position patients for lower-body examinations at a working height that is comfortable for the provider.



Elevated Flat Position
Easily position patients in a supine or Sims position at a working height that is more comfortable for the provider who prefers to work from a standing position.



BP Position
The Midmark® 627 with Patient Support Rails+ helps achieve proper patient positioning for blood pressure capture—arm supported at heart height, back supported and feet flat on the floor.

Better accessibility is better care.

If you're still using a standard, fixed-height examination table, chances are there have been patients you've had to examine while they sat in a side chair or wheelchair because they couldn't get onto the examination table. This compromises the exam—ensure that your patients receive equal access to quality care regardless of mobility. Ritter® examination chairs by Midmark adjust in height to make it easier for even the frailest patient to be seated.



Ritter 224
Barrier-Free®
Examination Chair



HELP PATIENTS HELP THEMSELVES.

A low seat height of 46 cm (18") and wide, flat transfer surface enable patients to transfer onto and off of the exam chair with minimal assistance.



STANDARD FEATURES THAT HELP ADD VERSATILITY



STIRRUPS + PELVIC TILT
Help position patients properly and comfortably for lower-body examinations.



CLEAN ASSIST ROLLER SYSTEM
The retractable roller base allows you to easily and safely move the chair for cleaning. (Available on the 627-011 model only.)



ACCESSORY RECEIVERS
Standard receivers are designed to make the integration of accessories easy and efficient—no tools required.



FOOT CONTROLS
Standard foot control designed to make patient positioning quick, effective and always within reach. (627-011 and 224-011 foot controls shown.)

ACCESSORIES

PATIENT SUPPORT RAILS+



Help achieve proper positioning of the patient's arm for blood pressure capture.

PATIENT SUPPORT RAILS



Provide patients with a 1¼ inch diameter continuous gripping surface for entering, exiting or repositioning on the exam chair.

ASSIST ARMS



Designed to offer patients a stable gripping surface for accessing the chair and a sense of security during positioning.

ARTICULATING KNEE CRUTCHES



The unique ball-socket design allows for a maximum range of adjustment, designed to ensure optimum knee and leg support. Available in a non-articulating option.

SUPPLY ASSISTANT



Provides extra storage for bulky exam supplies. Keeps drapes, gowns and similar items within easy reach.

STAINLESS STEEL TREATMENT PAN



Ideal for lower-body exams or other situations in which removal of fluid and debris is necessary.

UPHOLSTERY OPTIONS



SEAMLESS
Plush, seamless upholstery designed to be more comfortable for the patient and easier to clean.



FLAT
Provides an additional working surface for exams and procedures where you need a little more space.

PREMIUM THICK FOOTREST PAD



Provides additional support and comfort for patients in the flat position.

COLOR OPTIONS

STANDARD UPHOLSTERY COLORS

-  Mist 813
-  Stone 814
-  Restful Path 858
-  Iron Ore 851
-  Lunar Gray 845
-  Obsidian 857
-  Oasis 843
-  Soothing Blue 855
-  Dream 847
-  Mineral 844
-  Healing Waters 848
-  Shaded Garden 853
-  Citrus 849
-  Sandy Retreat 850
-  Curative Copper 846
-  Cranberry 859
-  Dark Linen 856
-  Latte 860
-  Robust Brown 852
-  Deep Earth 854

ULTRAFREE UPHOLSTERY COLORS

-  Stone 865
-  Restful Path 868
-  Obsidian 867
-  Harbor 861
-  Wheatgrass 862
-  Spice 863
-  Cranberry 869
-  Dark Linen 866
-  Latte 870
-  Branch 864

SPECIFICATIONS

Patient Weight Capacity:
500 lb (227 kg)

Height:
Uncompressed - 18" Low/37" High
(46 cm/94 cm)

Seamless Upholstery Dimensions:
28" W x 60" L
(71 cm x 152 cm)

Length with Footrest Extended:
76" (193 cm)

Pelvic tilt: 5 1/5°

Paper Roll Holder: Holds
24" x 3.56" dia. paper roll
(61 cm x 9 cm)

Chair Weight: 409 lb (185 kg)
Shipping Weight: 465 lb (211 kg)

Electrical Ratings:
224-011: 230 VAC, 50/60 Hz, 3.5 A
627-011: 230 VAC, 50/60 Hz, 5 A

ACCESSORIES:

- Patient Support Rails+
- Patient Support Rails
- Assist Arms
- Supply Assistant
- Articulating Knee Crutches
- Stainless Steel Treatment Pan
- Premium Thick Footrest Pad

627/224 COMPARISON CHART

	18" Seat Height (Uncompressed)	500 lb Weight Capacity	Power Base	Standard Rail Accessory Receivers	Active Sensing Technology*	Front Storage Drawer + Pelvic Tilt	Power Back	Clean Assist Roller System
627-011	•	•	•	•	•	•	•	•
224-011	•	•	•	•	•	•		



Designing better care.®



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Authorized CE Representative Countries in the EEC should direct all questions, incidents and complaints to the Midmark Authorized CE representative listed below:

CEpartner4U, Esdoornlann 13, 3951 DB Maarn, The Netherlands
Phone: +31.343.442.524 Website: cepartner4u.com

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

Use of TM and ® in this document denotes Midmark trademark rights in the United States. Midmark rights in these trademarks may also exist in other jurisdictions.

© 2021 Midmark Corporation, Miamisburg, Ohio USA
Products subject to improvement changes without notice. Litho in USA.

007-10072-00 Rev. D1 (9/21)

