



FAST. EFFECTIVE. COST-SAVING MAINTENANCE.

Quickly perform dental handpiece maintenance (cleaning and lubrication) for up to 4 handpieces at a time with the Lubricare 2[®]. Its clean air-blow system ensures thorough cleaning and lubrication inside your handpieces and consistent, repeatable high standards without the variability that may occur during manual lubrication.

FAST

Clean a turbine in approximately 20 seconds and a contra-angle handpiece in 40 seconds with preprogrammed operating modes that enable automatic maintenance.

EFFECTIVE

Achieve best cleaning and lubrication outcomes with separate dispensing of cleaner and lubricant for more precise application than systems that dispense the solutions together.

COST-SAVING

Reduce maintenance oil consumption by up to 4 times compared to manual maintenance procedures.

THE PERFECT PARTNER TO YOUR AUTOCLAVE

Using a handpiece oiler can help to save your practice money and time, as well as ensuring the maintenance of your handpieces is done correctly. The use of an automatic oiler such as the Lubricare 2 can stop excess oiling, reducing costs but also just as importantly protecting other equipment in practice.

The oil handpieces require if applied in excess can cause issues with autoclaves, the rubber components of your autoclaves such as seals and vacuum pumps can become brittle and perish leading to possible safety issues as well as expensive repair bills, the correct level of oil being used reduces the risk to your autoclave.

SPRAYNET[®] + LUBRIFLUID[®]

Use Spraynet and Lubrifluid for cleaning and lubrication to maintain impeccably hygienic instruments with an unrivaled service life.



Compatibility	E-type attachments (additional couplings may be required)
Cycle Duration	20 s turbine, 40 s contra angle
Capacity	4 handpiece lanes (2 cleaning, 2 lubricating) and 1 chuck lane
Rating	AC 240 V 50–60 Hz
Power Consumption	10 VA–25 VA
Air Flow Rate	60 NL/min–100 NL/min (0.3 MPa–0.5 MPa)
Air Pressure	43.5 psi–72.5 psi (0.3 MPa–0.5 MPa) [58 psi (0.4 MPa) recommended]
Weight Approx.	17 lb (8 kg)
Outer Dimensions	11.81" L x 11.81" W x 14.37" H (300 mm x 300 mm x 365 mm)

PLEASE VISIT OUR WEBSITE FOR PRODUCT INFORMATION



DENTALEX Ltd.

Singleton Court Business Centre
 Wonastow Road West
 Monmouth
 NP25 5JA

Phone: 01600 719 919

Email: info@dentalexmedical.co.uk

Web: www.dentalexmedical.co.uk

www.linkedin.com/company/dentalex-in-partnership-with-newmed

