



Tel: +86 755 2669 2869 Website: www.melevet.com E-mail: sales01@melevet.com

Shenzhen melevet Medical Co., Ltd.

ACCESSORIES



Y breathing circuit Coaxial breathing circuit



Waste gas hose



Oxygen hose



10





Endotracheal tube Non-rebreathing circuit

SD35 trolley

SD30 trolley









Anesthesia face masks

Breathing Bags

Endotrachealtube support Oxygen face masks





AM70

VETERINARY ANESTHESIA MACHINE

- Ultra-compact design & multiple installation methods save more space of operating rooms.
- 10 inches screen, clear observation. touch screen, smooth and simple operation.
- Minimum 5ml tidal volume better supports small animals.
- Turbine drive, no need to configure drive gas, realize low-cost anesthesia.
- Integrates 2 auxiliary power outputs for easy device integration.
- Support VCV, PCV, VS, SIMV multiple modes.

AM30

VETERINARY ANESTHESIA MACHINE

- Ultra-compact design & multiple installation methods save more space of operating rooms.
- Selectatec vaporizer, quick change.
- ACGO switch, switch between rebreathing anesthesia and non-rebreathing circuits in one step.
- Equipped with oxygen pressure gauge and safety valve to ensure the safety
- Pre-calibrated APL valve precisely sets maximum airway pressure.
- 1400ml large capacity CO2 absorber canister.





AM10

VETERINARY ANESTHESIA MACHINE

- Ultra-compact design & multiple installation methods save more space of operating rooms.
- Oxygen pressure gauge and safety valve ensure the safety.
- Selectatec vaporizer, quick change.
- Up to six channels can be flexibly configured.

VR30

SERIES VAPORIZERS

- Selectatec and Cagemount two installation methods.
- Sevoflurane, Enflurane, Isoflurane optional.
- Keyfill, Easyfill, Quickfill, Pourfill optional.
- Flow, pressure and temperature compensation ensure accurate output.



03



BN30

BAIN BLOCK

- High precision airway pressure gauge, 0-100cmH20.
- Pre-calibrated APL valve to accurately set airway pressure in one step.
- The APL button can quickly increase the airway pressure by about 30cmH20.
- Perfect fit for AM30 anesthesia machine.

PA50

PATIENT PRESSURE ALARM

- Multiple alarm gears can be set with an accuracy of 5%.
- Standard 9V battery.
- Small and compact size, flexible and convenient installation.
- Low battery alarm.





PA30

OXYGEN PRESSURE ALARM

- Automatic alarm below 200kPa
- Standard 9V battery.
- Small and compact size, flexible and convenient installation.
- Low battery alarm.

05